

Over-expression and Rescue Constructs (Yuan Lab)



Note 1: This protocol is based on the pEarleyGate vectors (ABRC Stock CD3-683 – CD3-695, CD3-724)

1. PCR insert and purify

use a high-fidelity enzyme (Phusion) to amplify your insert

<u>PCR</u>	<u>1x</u>	
dH ₂ O	12.0 μ l	
5x Phusion buffer	4.0 μ l	
dNTPs (10mM)	0.4 μ l	
DMSO 100%	0.6 μ l	
F primer (5 μ M)	1.0 μ l	*add "CACCC" at beginning of forward primer
R primer (5 μ M)	1.0 μ l	
template	1.0 μ l	
<u>Phusion enzyme</u>	<u>0.2 μl</u>	<u>ADD THIS LAST</u>
TOTAL	20.0 μ l	

PCR Program:

1. 98°C for 30 sec
 2. 98°C for 10 sec
 3. 60°C for 15 sec
 4. 72°C for 30 sec-2min
 5. go to step #2, 33x
 6. 72°C for 10 min
 7. 12°C forever
- *anneal temperature depends on the primers
*depends on the length of the fragment (30 sec/1kb)

Run PCR product on a 1% gel to check for size and a single band. If the PCR product is a clean, single band, it can be PCR purified. If there are multiple bands, gel purify the fragment. Elute PCR product in 25 μ l dH₂O.

2. pENTR-D TOPO cloning reaction

dH ₂ O	1.0 μ l	
PCR product	1.0 μ l	*adjust the volume of PCR product and dH ₂ O based on concentration.
salt solution	0.5 μ l	
Topo vector	<u>0.5 μl</u>	
	3.0 μ l	

Mix and incubate at room temperature for at least 30 minutes, then put on ice.
Transform entire reaction into *E.coli* competent cells and plate all on Kan plates.

Transforming One Shot Top10 *E. coli* competent cells:

- a. Take out S.O.C. from -20 freezer to thaw
- b. Thaw 1 vial of One Shot cells on ice (1 vial can do 3 or 4 transformations, but it cannot be re-frozen) and aliquot the cells into separate PCR tubes, if necessary
- c. Add product to be transformed into One Shot cells and mix gently (DO NOT mix by pipetting up and down)
- d. Incubate on ice for 5 minutes
- e. Heat shock the cells for 30 seconds at 42° C in the PCR machine
- f. Immediately transfer the tubes to ice
- g. Add 250 µl S.O.C. to the tube (125 µl if the cells were split into 2 or 3 tubes)
- h. Shake the tube horizontally at 37° C for 1 hour
- i. Spread 50-200 µl of the transformation on pre-warmed selective plates and incubate overnight at 37° C

3. Colony PCR to check for the insert

Use **M13F (or T7F)** primer and an insert specific primer (**cdsR**) to check for the size and direction of the insert – **Important!** Sometimes the insert could be inserted in the opposite direction

Making a replica plate and adding template:

- a. Circle and number the colonies you wish to colony PCR (8-16 colonies)
- b. Get a new selective plate and make a grid and number each cell = replica plate
- c. Using a P20 pipet tip, gently touch the side of a SINGLE colony, gently touch the corresponding grid on the replica plate, then place the pipet tip in the PCR tube containing your Master Mix
- d. Incubate replica plate at 37° C for a few hours or overnight

Run PCR product on a gel to check for the presence of an insert in the correct orientation.

4. Culture two colonies with the correct insert and isolate the plasmid = entry clone

In a labeled 15mL Falcon tube, pipet 3 mL of LB broth.

Add 3 uL of Kan to each tube

Pick a colony with the insert with a P20 tip and eject the tip into the Falcon tube.

Incubate at 37°C with shaking overnight (try to start your cultures in the afternoon and take them out of the shaker in the morning – if left to grow too long, it will start to die).

Isolate the plasmid from the culture using a Plasmid Mini-prep Kit. Elute plasmid in 70 uL of dH₂O. **note: elute the plasmid twice (35 + 35) works much better than once (70).**

5. PCR using M13 F/R primers to amplify fragment for LR reaction

Necessary since the entry clone and the destination vector both have Kan resistance.
Use the Phusion protocol (from Step 1) and M13F/R primers to amplify the fragment.
Use only 28 cycles since the PCR should be very strong
Gel or PCR purify the fragment.

6. LR Recombination reaction

dH ₂ O	2.5 ul	
entry clone PCR fragment	1.0 ul	
Destination vector	0.5 ul	*select desired vector p100, p103, p302, etc.
LR clonase	<u>1.0 ul</u>	*vortex 2 sec. twice and spin down before adding
	5.0 ul	

Mix well and incubate at room temperature (25°C) for 2 hours.
Add 1 ul proteinase K to terminate the reaction and vortex briefly.
Incubate at 37° for 10 minutes.
Transform 3 ul of the reaction into *E. coli* competent cells and plate 20-100 uL on Kan plates.

7. Colony PCR using insert specific primers to check for insert

Use insert specific primers to check for insert or a primer on the vector if available (attR2 for example)
Make a replica plate.
Run PCR product on a gel to check for the presence of an insert.

8. Culture two colonies and isolate the plasmid = final plasmid

Culture two colonies in 3 mL of LB broth + Kan overnight at 37°C with shaking.
Isolate the plasmid using a Plasmid Mini-prep Kit. Elute in 70 uL (35 + 35) of dH₂O.

9. Sequence verify the final plasmid