

RNAi Plasmid Construction Using pFGC5941 (Yuan Lab)

Note 1: This protocol is based on the vector pFGC5941 (ABRC Stock CD3-447).

Note 2: To avoid off-target effect, make sure no other regions in the interested genome perfectly match the RNAi fragment (150-500 bp) for a contiguous block longer than 16 bp. Also, make sure there are no restriction sites for the enzymes NcoI, AscI, BamHI, or XbaI within the RNAi fragment.

Note 3: When designing primers to amplify the RNAi fragment. Add “GTTCTAGACCATGG” at the 5’ end of the Forward primer and add “GTGGATCCGGCGCGCC” at the 5’ end of the Reverse primer.

Note 4: Make sure you have digested the pFGC5941 vector using NcoI/AscI before the first ligation.

Protocol:

1. Amplify insert from cDNA or gDNA (if the fragment contains no intron) using Phusion PCR;

Do **TWO** 20- μ l reactions:

4 μ l 5x Phusion buffer

0.5 μ l 10mM dNTPs

0.6 μ l DMSO

1.0 μ l template

0.2 μ l Phusion enzyme

11.0 μ l dH₂O

1.5 μ l 5 μ M primer F

1.5 μ l 5 μ M primer R

20 μ l total

Phusion PCR program:

cycle 1: 98°C for 0:30

cycle 2: (32x) 98°C for 0:10

58°C for 0:20 (or whatever ideal annealing temperature)

72°C for 0:30

cycle 3: 72°C for 5:00

cycle 4: 12°C for ever

2. Digest insert with NcoI/AscI and BamHI/XbaI

2.5 μ l 10x CutSmart Buffer

4.5 μ l dH₂O

1.5 μ l NcoI

1.5 μ l AscI

15 μ l PCR product

25 μ l total

same protocol for BamHI/XbaI digestion;

incubate 37°C for 1 hour;

gel purify digests and save the BamHI/XbaI digested insert for the second ligation.

3. First ligation (Want an insert to vector molar ratio of 2:1 to 6:1)

2 μ l linearized pFGC5941 digested with AscI/NcoI (~175ng; **adjust volume as needed**)
4 μ l insert digested with AscI/NcoI (~15-30ng)
2 μ l T4 ligase buffer
1 μ l T4 ligase
11 μ l dH₂O
20 μ l total

incubate 30 minutes at room temperature;
transform 10ul into *E. coli* competent cells (homemade) and plate on Kan plates.

4. Colony PCR to check for first insert

Circle the biggest colonies on your plate and label them 1-8.

Make a replica plate for your colonies.

PCR across the first insert using primers on the vector to check for an insert:

An empty vector will give a band of 700bp

8.0 μ l dH₂O
1.0 μ l 10x buffer
.125 μ l dNTPs
0.5 μ l pFGC5941 **2372 F**
0.5 μ l pFGC5941 **3082 R**
0.05 μ l Taq
10 μ l total

Colony PCR Program:
cycle 1: 95°C for 3:00
cycle 2: (32x) 95°C for 0:15
55°C for 0:15
72°C for 1:00
cycle 3: 72°C for 7:00
cycle 4: 12°C forever

5. Pick two correct colonies and inoculate into 3 mL LB+Kan broth

incubate in 37°C shaker overnight

The next day, do a plasmid prep (mini-prep kit) with 1 of the colonies that grew well

6. Digest plasmid with BamHI/XbaI

5 μ l 10x CutSmart Buffer
12 μ l dH₂O
1.5 μ l XbaI
1.5 μ l BamHI
30 μ l plasmid * adjust volume based on concentration; you want 2000-5000 ng of plasmid
50 μ l total

37°C for 1 hour
gel purify digest

7. Ligation #2

2 μ l vector that contains the first insert, digested with BamHI/XbaI (~175ng; **adjust volume based on concentration**)

4 μ l insert digested with BamHI/XbaI (done in step 2) (want ~15-30 ng)
2 μ l T4 ligase buffer
1 μ l T4 ligase
11 μ l dH₂O
20 μ l total

incubate 30 minutes at room temperature
Transform 10ul into *E. coli* competent cells (homemade) and plate on Kan plates

8. Colony PCR to check for second insert

pFGC5941 **3930 F** & pFGC5941 **4430 R**

Vector without insert will give a band of 500bp

9. Pick two correct colonies and inoculate into 3 mL LB+Kan broth

incubate in 37 degree shaker overnight

Plasmid prep (mini-prep kit)

10. Check plasmid for inserts

PCR to check for both inserts:

2372F/3082R or RNAi_R (insert specific)

3930F/4430R or RNAi_F (insert specific)

11. Sequence to verify

Use 4 primers:

2372F, 3082R, 3930F, 4430R

Note: in the sequencing reaction, add DMSO to aid in the sequencing across the restriction enzyme digest sites (the chromatogram peaks usually drop off dramatically right after the digest sites; an alternative strategy is to PCR the final plasmid with 2372F&3082R for the left insert and 3930F&4430R for the right insert and then sequence the PCR product)

12. Transform into agrobacterium for infiltration

Primer sequences:

pFGC5941_2372F: CTTCATCGAAAGGACAGTAGAA

pFGC5941_3082R: CCAAACAGGCTCATAGATACT

pFGC5941_3930F: TGTACATCAGAATGTTTCTGAC

pFGC5941_4430R: CGCTCTATCATAGATGTCGCTA